

# Update

Dear Practice Nurses and Doctors,

We have changed our collection requirements for Zinc, Copper and other trace elements.

## As of 27<sup>th</sup> July 2023:

For Zinc and Copper testing we are now utilising  
a new Trace Elements tube  
(Dark Blue with Red stripe on label)  
which contains a clot activator coating.



For all other trace elements,  
we are utilising EDTA tube (Lavender).



## Doctors Order Request

Surgery :	
Address	
Contact Name:	Phone:
Date:	Received: (office use)

### Doctor's request for stores GOODS TO BE USED FOR 4CYTE PATHOLOGY COLLECTIONS ONLY.

Under Medicare Australia guidelines, collection equipment is only supplied if used for pathology collection. Fortnightly ordering is recommended for efficient stock rotation. Please allow 5 working days for delivery.

DESCRIPTION (UNIT OF ISSUE)	QTY	MAX	DESCRIPTION (UNIT OF ISSUE)	QTY	MAX
<b>VACCUUM BLOOD COLLECTION</b>			<b>URINE/STERILE COLLECTION</b>		
Barrel for multisample needles (Bag of 250)		1	Specimen Jar (Yellow top) (Each)		50
Needle: 21G X 1 1/2" multisample (Box 100)		1	24hr Plain Urine Collection Container: Nil preservative (Each)		2
Needle: 22G x 1 1/2" multisample (Box 100)		1	24hr Acid Urine Collection Container (Each)		2
Needle: 21G Butterfly Green (Each)		10	24hr Urine Collection Funnel (Each)		2
Needle: 23G Butterfly Blue (Each)		10	Paediatric Urine Bag (Each)		10
Tube: EDTA (Lavender top) 4ml (Each)		100	<b>GYNAECOLOGICAL AND CYTOLOGICAL COLLECTION</b>		
Tube: SST (Gold top) 8.5ml (Each)		100	ThinPrep® (Tray 25)		2
Tube: Sodium Citrate (Light Blue top) 2.7ml (Each)		100	Rovers/Wallach Cervex Sampler (Pack 25)		2
Tube: EDTA (Blue top) Trace Element 6ml (Each)		6	Endo-Cervical Brush (Pack 10)		2
Tube: Fluoride Oxalate (Grey top) 2.0ml (Each)		100	FLOQ Swab Self Collect CST (Each)		10
Tube: Lithium Heparin (Green top) 4ml (Each)		6	Speculum – EOS Small (Each)		20
Tube: EDTA Gel PPT (White/Pearl top) 5ml (Each)		6	Speculum – EOS Medium (Each)		20
Tube: Plain Clot (Red top) 6ml (Each)		6	Speculum – EOS Large (Each)		20
Microtainer Tube: EDTA (Lavender) 1ml (Each)		1	Speculum – VIKING Small (Each)		20
Microtainer Tube: SST (Gold) 1ml (Each)		1	Speculum – VIKING Medium (Each)		20
Microtainer Tube: Fluoride Oxalate (Grey) 1ml (Each)		1	Speculum – VIKING Large (Each)		20
Tube: Citrate (Light Blue) 1.8ml (Each)		1	Cotton Applicator: Single tip: 15cm (Pack 100)		1
Blood Culture Bottle – Aerobic (Blue Cap) (Each)		1	<b>HISTOPATHOLOGY</b>		
Blood Culture Bottle – Anaerobic (Purple Cap) (Each)		1	Histology Container: 25ml with formalin (Tray 24)		1
Blood Culture Bottle – Paeds (Pink Cap) (Each)		1	Histology Container: 70ml with formalin (Tray 10)		1
Kit: Percept NIPT (VCGS)		4	Biopsiblade® (Box 25)		1
Kit: Quantiferon TB Gold Plus (Set of 4 tubes)		5	Biopsy Punch: 3mm with plunger (Box 20)		1
<b>NON-EVACUATED BLOOD COLLECTION</b>			Biopsy Punch: 4mm with plunger (Box 20)		1
Syringe: 5ml without needle (Box 100)		1	Biopsy Punch: 6mm (Box 20)		1
Syringe: 10ml without needle (Box 100)		1	Biopsy Punch: 8mm (Box 20)		1
Needle: 21G x 1 1/2" (Box 100)		1	<b>STATIONERY</b>		
Needle: 23G x 1 1/4" (Box 100)		1	A4 Pathology Request Form with labels (Pack 100)		5
<b>FAECES COLLECTION</b>			Pathology Request Pad (Carbon) (Book)		2
Faeces Jar (Brown top) (Each)		20	A4 Maternal Serum Screen Form (Pad) VCGS		1
Kit: FOB (Set of 3 Jars)		10	A4 Percept NIPT Form (Pad) VCGS		1
<b>SWAB COLLECTION</b>			A4 Prepair Genetic Carrier Screening Form (Pad) VCGS		1
eSwab Regular Collection (Blue cap) (Each)		20	Doctor's Request for Stores (Sheet) (Each)		10
eSwab Minitip (NPG/Urethral) (Orange Cap) (Each)		20	Swab and Brush Guide		6
<b>GENERAL</b>			Tube Guide		6
Clear Specimen Bag (Pack 100)		5	<b>MISCELLANEOUS</b>		
Red Specimen Bag (Pack 100)		2			
Glucose Load Drink 75gm (Each)		5			
Alcohol Swab (Pack 100)		3			
Urea Breath Test Kits C14 (Each)		5			